



Avoid misfiring when adjusting the depth. Ergonomic design for improved holding.

Improved safety cap removal – cap automatically semi ejects when turned to remove.



SAFETY LANCETS

Vitrex® Flex 3 Depth Adjustable Safety lancet

The All-In-One for emergency use.

USER CONVENIENCE adjustable design provides three different penetration depths from a single unit - easy to use pre-loaded and contact activated device. Trigger position indicates the before and after use.

USER SAFETY is assured because the needle is fully shielded before and after use.

COMFORT is provided by high quality tri-bevel needle with high-speed penetration.

PATIENT SAFETY sterile, single-use, unable to be re-used.

APPLICATION can be used for blood glucose and other diagnostic tests.


Technical Specs.

Colour Code	Specification	Article No.
Blue	23G 1.3 mm - 1.8 mm - 2.3 mm Penetration Depth	395813
Yellow	26G 1.3 mm - 1.8 mm - 2.3 mm Penetration Depth	396813
Orange	28G 1.3 mm - 1.8 mm - 2.3 mm Penetration Depth	397813

Shelf life 5 years



Packaging:
Primary packaging: 100 pcs.
Transit packaging: 1,000 pcs.
Pallet: 84,000 pcs.



Penetration depth adjustment is done without touching the activation bottom.